

## **Briteline**

## **Product Data Sheet**

**Product Description: Briteline Urine Temperature Indicators** 

## Standard Range: 89 to 101°F or 31.5°C to 38.5°C

## **Construction:**

Adhesive Tape	Adhesive: Single-sided adhesive paper label (Fasson S- 246)
Covering Film	Polyester
Color Change Material	Non-toxic mixture of liquid crystal esters
Accuracy	L215F ±1°F L215C ±1°C L216F ±1°F
Warranty Period	12 months from invoice date when stored in optimal storage conditions (between 60 and 80°F and 40-60% relative humidity)

Application Data: Apply temperature indicating label to outside of urine specimen cup. Reversible analog line moves along the scale to show detailed temperature reading in °F or °C.
Function: Display of temperature, as indicated by a bright line, of urine specimen cup. Label temperatures are adjusted with a 1°F / 0.6 °C offset.
Use: Ideal for use in clinics, hospitals, and doctors' offices.