

## **Clinitemp™ Product Data Sheet**

**Product Description:** L314 Latex Free Forehead Moving Line Temperature Monitor

**Standard Range:** 84-106°F / 29-41°C

**Construction:** 

Adhesive Tape Adhesive: Acrylic copolymer pressure sensitive adhesive, thickness 37µ each side of carrier

Carrier: Non-woven tissue, thickness 33µ Release Liner: Fully bleached, thickness 178µ Operating temp: -1 to 82°C (continuous use)

Covering Film Polyester

Color Change Material Nontoxic mixture of liquid crystal esters

Accuracy  $\pm 1.00^{\circ} \text{F} / 0.6^{\circ} \text{C}$ 

Shelf Life 12 months from invoice date when stored at room temperature & humidity

(i.e. 21°C (70°F) & 50% relative humidity).

**Application Data:** Peel temperature)indicating label from backing paper. Apply to forehead. Reversible analog

line moves along the scale to show detailed temperature readings in °F and °C.

**Function:** Continuous display of temperature, as indicated by a bright line, of anesthetized patients.

Core adjusted with 3.5°F / 2°C offset.

**Use:** Ideal for use in Pre-op, Surgery, PACU

Hallcrest

www.spotsee.io Rev: 03/2021